




ASSESSING REAL-TIME SICKLING KINETICS IN INDIVIDUAL CELLS FROM ERYTHROCYTE POPULATIONS OF SUBJECTS WITH DIFFERENT HbA, HbS AND HbF PERCENTAGES

ALIYA U. ZAIDI, PHD¹, XIUFENG GAO, MD, MSC¹, ROBERT GOODRICH, BSC¹, MARTA FERRANTI, MSC¹, RASA BORHAN, BSC¹, KE LIU, PHD^{1*}, NNAMDI OKEKE, PHARMD¹, MICHAEL TARASEV, PHD, MBA¹, MARIAMA KABORE², COURTNEY D. FITZHUGH, MD², PRIYA S. CHOCKALINGAM, PHD³ AND PATRICK C. HINES, MD, PHD^{1,4}.

¹Functional Fluidics Inc., Detroit, MI ²NHLBI/NIH, Bethesda, MD ³Beam Therapeutics, Cambridge, MA ⁴Wayne State University Division of Pediatric Critical Care Medicine, Detroit, MI



INTRODUCTION

- Polymerization of deoxygenated hemoglobin S (HbS) leads to red blood cell (RBC) sickling and vaso-occlusion
- Allogeneic hematopoietic stem cell transplantation (HSCT) is the only approved curative option available for sickle cell disease (SCD)
- Objective biomarkers are needed to assess disease status, population-scale sickling kinetics and to monitor the pharmacodynamic impact of Hb-modifying and curative therapies
- Functional Fluidics has developed a dynamic sickling assay (DSA) that uses an enzymatic O₂ scavenging system providing tight control over the rate and depth of induced-hypoxia

AIM

Assess real-time RBC sickling kinetics in sickle cell trait (SCT) and SCD patients with variable distributions of HbA, HbS and HbF, and compare pre-and post-HSCT

METHOD

Dynamic Sickling Assay (Figure 1)

- Hypoxia is induced through a timed protocatechuic acid/protocatechuate-3,4-dioxygenase enzymatic reaction (PCD-PCA)
- The fraction of sickled cells is determined by time-lapse imaging, analyzed with AI technology and reported as a function of time
- Measured parameters include:
 - **mPoS@5%**: time to reach 5% induced sickling (min)
 - **mPoS@50%**: time to reach 50% induced sickling (min)
 - **Rate of Sickling**: highest induced sickling rate (%·minute⁻¹)
 - **Maximum Sickling**: highest sickling percentage (%)
 - **AUC at 10min**: The area under the curve during 10min of DSA (%·minute)

Sample collection

- SCT and SCD subjects, with and without chronic exchange transfusions were enrolled in The Institute for Regenerative and Cellular Medicine (IRCM) IRB under the RBC health initiative protocol (FF-RBC-003v5)
- SCD subjects were enrolled in NIH haploidentical HSCT protocol 17-H-0069. Blood samples from 4/10 planned subjects were collected pre-and 3, 6, and 12 months after HSCT




Figure 1. Dynamic Sickling Assay setup

RESULTS

SICKLING KINETICS IMPROVE AS HbS% DECREASES AND HbF% INCREASES




Figure 2. Sickling kinetics in SCT subjects (A), transfused SCD subjects (B) and SCD subjects with variable HbF (C)

CONCLUSIONS

- SCT subjects with similar HbS% exhibit different sickling kinetics upon exposure to the same hypoxic stress due to intracellular RBC distribution of HbA and HbS
- Transfusion SCD samples show an improvement in sickling kinetics as HbS% decreases approaching a significant correlation with AUC10
- Statistically significant correlations are observed with increasing HbF percentages on 4 out of 5 sickling parameters
- Compared to pre-transplant values, post-HSCT samples exhibit significant improvement in: mPOS@5%, rate of sickling, maximum sickling and AUC10
- DSA can differentiate samples based on Hb distribution and can be used to assess sickling kinetics pre- and post HbS modifying and curative therapies, such as hematopoietic stem cell (HSCT)

RESULTS

POST-HSCT SAMPLES EXHIBIT IMPROVED SICKLING KINETICS




Figure 3. Sickling kinetics in SCD subjects pre-and post HSCT (A) Representative sickling profile from 1 patient (B) Comparative analyses from 4 SCD subjects *only 1 post-HSCT sample exhibits >50% max. sickling

- Compared to baseline values, post-HSCT samples exhibit statistically significant improvement in: mPOS@5%, rate of sickling, maximum sickling and AUC10

Table 1. Sickling kinetics correlations

		mPoS@5%		mPoS@50%		Rate of Sickling		Max Sickling		AUC10	
	n	r	p	r	p	r	p	r	p	r	p
HbS% (Sickle Cell Trait)	3	NA	NA	NA	NA	NA	NA	0.30	0.803	0.86	0.345
HbS% (Transfused SCD)	3	-0.99	0.100	NA	NA	0.626	0.596	0.99	0.107	0.99	0.066
HbF% (non-Transfused SCD)	9	0.87	0.002	0.83	0.006	-0.74	0.023	0.15	0.692	0.73	0.025

- No significant correlation in SCT subjects, although noteworthy differences evident under hypoxic conditions in samples with similar HbS% levels
- HbS% and DSA parameters show a correlation in SCD transfused samples, not statistically significant likely due to limited sample size
- HbF% demonstrates a significantly strong correlation with DSA parameters except for max. sickling

REFERENCES

- Galal, A., & Asslan, M. (2019). Allogeneic hematopoietic stem cell transplantation in sickle cell disease. *Current Opinion in Hematology*, 26(6), 399-405.
- Tarasev, M., Gao, X., Ferranti, M., Goodrich, R., Edenstrom, A., Zaidi, A., Hines, PC. Dynamic Sickling Assay: Using Enzymatically Induced Hypoxia to Assess Real-Time Sickling Scalable from Large Erythrocyte Populations to Individual Cells. *Blood*, 2022. 140 (Supplement 1): p. 4976-4977.
- Kato, G.J., Piel, F.B., Reid, C.D., Gaston, M.H., Ohene-Frempong, K., Krishnamurti, L., Smith, W.R., Panepinto, J.A., Weatherall, D.J., Costa, F.F. and Vichinsky, E.P., 2018. Sickle cell disease. *Nature Reviews Disease Primers*, 4(1), pp.1-22.

ACKNOWLEDGEMENTS



CONTACT INFORMATION

